

	Test	Uncorrected J1 or better (20/25)	Uncorrected J3 or better (20/40)	Source
MA60BM* (Monocular)	Near Card (33cm)	0.6%	11.9%	AcrySof ReSTOR Package Insert
CeeOn91 I A^ (Monocular)	Near Card (33 cm)	0.9%	6.2%	Tecnis MF Package Insert
Softec HD (Monocular)	Rosenbaum Pocket Screener (14")	4.1%	27.9%	James P. Gills, MD <i>ASCRS 2011</i>

	Test	Uncorrected J1 or better (20/25)	Uncorrected J3 or better (20/40)	Source
MA60 BM* (Monocular)	Near Card (33cm)	0.6%	11.9%	AcrySof ReSTOR Package Insert
CeeOn91 I A^ (Monocular)	Near Card (33 cm)	0.9%	6.2%	Tecnis MF Package Insert
Softec HD (Monocular)	Rosenbaum Pocket Screener (14")	5.9%	35.3%	T. Hunter Newsom, MD <i>ASCES SEE Meeting 2011</i>

	Test	Uncorrected J1 or better (20/25)	Uncorrected J3 or better (20/40)	Source
MA60BM* (Binocular)	Near Card (33cm)	2.5%	26.1%	AcrySof ReSTOR Package Insert
CeeOn91 I A^ (Binocular)	Near Card (33 cm)	1.7%	21.0%	Tecnis MF Package Insert
Softec HD (Binocular)	Rosenbaum Pocket Screener	37.9%	57.6%	J. Pit Gills, MD <i>ASCES SEE Meeting 2011</i>

		Uncorrected J1 or better (20/25)	Uncorrected J3 or better (20/40)	
	Test			Source
MA60BM* (Binocular)	Near Card (33cm)	2.5%	26.1%	AcrySof ReSTOR Package Insert
CeeOn91 I A^ (Binocular)	Near Card (33 cm)	1.7%	21.0%	Tecnis MF Package Insert
Softec HD (Binocular)	Rosenbaum Pocket Screener	37.9%	57.6%	J. Pit Gills, MD <i>ASCES SEE Meeting 2011</i>